

Check the product label for actual catalog number, lot and expiry date.

## Taq DNA Polymerase

CAT.#	SIZE	COMPONENTS	COMPONENT COMPOSITION
PCE0201	1500 u	1500 u - Taq DNA Polymerase, 5 u/μl 4 x 2 ml - 10X PCR Buffer 2 x 2 ml - 50 mM MgCl <sub>2</sub>	Enzyme in storage buffer. 10X PCR Buffer contains enhancers and stabilizers, but no dNTPs and no Mg <sup>2+</sup> .
PCE0202	3000 u	2 x 1500 u - Taq DNA Polymerase, 5 u/μl 8 x 2 ml - 10X PCR Buffer 4 x 2 ml - 50 mM MgCl <sub>2</sub>	Enzyme in storage buffer. 10X PCR Buffer contains enhancers and stabilizers, but no dNTPs and no Mg <sup>2+</sup> .

Storage In the dark at -20°C.

### APPLICATIONS

- Routine PCR up to 5 kb
- RT-PCR
- Colony PCR
- TA cloning, library construction
- Genotyping, screening

### PRODUCT DETAILS

highQu Taq DNA Polymerase is the classical enzyme for routine PCR applications providing high amplification yields of 3-5 kb targets under various conditions.

Taq DNA Polymerase is purified from a recombinant *E. coli* strain carrying the *Taq* DNA polymerase gene.

Taq DNA polymerase is thermostable 5' → 3' DNA polymerase.

It lacks 3' → 5' exonuclease (proofreading) activity and has low 5' → 3' exonuclease activity. Polymerase exhibits

deoxynucleotidyl transferase activity resulting in A-overhang at the 3'-ends of PCR products and allowing for TA cloning.

The PCR accuracy of Taq DNA Polymerase is 4.5 x 10<sup>4</sup>

(nucleotides incorporated before the error occurs) and it produces A-tailed products suitable for ligating into TA cloning vectors.

### BENEFITS

- High yields in routine PCR and under fast cycling conditions
- Guaranteed successful DNA preparation for TA cloning
- Robust on complex templates

### MAGNESIUM AND NUCLEOTIDE USE

The 50 mM MgCl<sub>2</sub> is provided in a separate tube what allows for magnesium optimization. The dNTPs in mixes or sets can be purchased separately.

- Typical concentration of each dNTP in the reaction is 0.2 – 0.25 mM. Higher concentration increase yields, however Mg<sup>2+</sup>-ions bind to dNTPs, therefore, both components shall be present in coordinated concentrations. Too high dNTPs and magnesium concentrations reduce PCR fidelity.
- Mix well dNTP and magnesium solutions, to avoid concentration fluctuations. Use final 3 mM MgCl<sub>2</sub> and 0.25 mM each dNTP concentrations for routine PCR.

Starting dNTP Mix conc.	Vol. of dNTP mix per 50 μl r.	Final Mg <sup>2+</sup> conc. in r.	Vol. of 50 mM MgCl <sub>2</sub> per 50 μl rxn to achieve desired conc.
10 mM	1.25 μl	2 mM	2 μl
25 mM	0.5 μl	3 mM	3 μl

### PROTOCOL

- Take typical measures to prevent PCR cross over contamination, keep your bench clean, wear gloves, use sterile tubes and filter pipet tips.
- Include a no-template control and positive control in parallel.
- Thaw and keep reagents on ice. Mix well before use.
- The longer the amplicon, the longer the extension time. Taq DNA Polymerase speed is ~ 2000 nucleotides/min, so 15-90 seconds of extension can be provided per cycle, depending on amplicon size.
- Start annealing with 55°C and perform gradient by increasing temperature in 2°C up to 65°C to choose the best. Calculate primer annealing temperature using software.

✓ Prepare a 50 μl reaction:

Rev. & For. Primers	0.1-0.5 μM final (1-2 μl of 10 μM each)
cDNA Template <i>or</i>	<100 ng <i>or</i>
gDNA Template	5-500 ng
10X PCR Buffer	5 μl
dNTP (10 mM dNTP Mix NUM0201)	0.25 mM final (1.25 μl of 10 mM dNTP mix or 0.5 μl of 25 mM dNTP mix)
50 mM Mg Cl <sub>2</sub>	3 μl
Water (PCR Water WAT0110)	to 49 μl
Taq DNA Polymerase, 5 u/μl	0.25- 1 μl

✓ Mix gently, avoid bubbles. Place into the instrument set:

Initial denaturation	1 cycle: 95°C – 60 sec
Denaturation	30-40 cycles: 95°C - 15 sec
Annealing	30-40 cycles: 55-65°C – 15 sec
Extension	30-40 cycles: 72°C – 15-90 sec
Final extension	1 cycle: 72°C – 5 min (for TA cloning)

✓ Store probes for short time on ice, for long at -20°C.

IN VITRO RESEARCH USE ONLY

#### ORDERING

T: +49 7250 33 13 401  
F: +49 7250 33 11 413  
[order@highQu.com](mailto:order@highQu.com)  
[www.highQu.com](http://www.highQu.com)

#### SALES

T: +49 7250 33 13 401  
F: +49 7250 33 11 413  
[sales@highQu.com](mailto:sales@highQu.com)

#### TECHNICAL SUPPORT

T: +49 7250 33 13 401  
F: +49 7250 33 11 413  
[tech-support@highQu.com](mailto:tech-support@highQu.com)